Dilution Instructions for Wescor Concentrated Hematology Reagents for use with the Aerospray® Model 712#
English

Required Supplies and/or Equipment:

SS-135S: Surfactant (14 mL)
AC-043: 500 mL bottle
Graduated cylinders
Filtered Deionized Water
Anhydrous, Reagent Grade Methanol
Anhydrous, Reagent Grade Ethanol or Reagent Alcohol

If bottles are reused be sure to clean the bottle thoroughly. Use 1/10 mixture of household bleach to disinfect the bottles.

WARNING! These reagents contain moderately toxic chemicals that require care in handling. Always use appropriate safety measures including gloves and eye protection when handling reagents. Refer to material safety data sheet (MSDS).

SS-135/149A Hematology Reagent A Concentrate (135 mL)

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-035A):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL of SS-135/149A to 500 mL deionized water.
3. Add 1.5 mL SS-135S.
4. Mix well.
5. Label new bottle with supplied label.
6. Write LOT number and expiration date from concentrate SS-135/149A bottle on new label.

SS-135/149B Hematology Reagent B Concentrate (135 mL)

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-035/049B):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL of SS-135/149B to 500 mL deionized water.
3. Add 1.5 mL SS-135S.
4. Mix well.
5. Label new bottle with supplied label.
6. Write LOT number and expiration date from concentrate SS-135/149B bottle on new label.

SS-135/149C Hematology Reagent C Concentrate (135 mL)

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-035C):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL of SS-135/149C to 500 mL deionized water.
3. Add 1.5 mL SS-135S.
4. Mix well.
5. Label new bottle with supplied label.
6. Write LOT number and expiration date from concentrate SS-135/149C bottle on new label.

SS-135/149A Hematology Reagent A Concentrate (135 mL) prepared as Basofix™ program.

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-049A):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL of SS-135/149A to 475 mL deionized water, 12.5 mL ethanol, 12.5 mL methanol.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate SS-135/149A bottle on new label.

SS-135/149C Hematology Reagent C Concentrate (135 mL) prepared as Basofix™ program.

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-049C):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL of SS-135/149C to 425 mL deionized water, 37.5 mL ethanol, 37.5 mL methanol.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate SS-135/149C bottle on new label.

SS-149P Basofix™ Predip Concentrate (15 mL)

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-049P):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL of SS-149P to 500 mL methanol.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate SS-149P bottle on new label.

SS-148 Aerofix® Additive for Methanol (135 mL)

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-048):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL SS-148 to 500 mL Methanol.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-029C Nozzle Cleaning Solution Concentrate (250 mL)

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-029):

1. Add 250 mL Methanol or Ethanol to bottle of SS-029C.
2. Mix well.
3. Check box on label when mixed.

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