Wescor Gram Reagent Concentrate Dilution Instructions English

Required Supplies and/or Equipment:
AC-043 – empty 500 mL bottle
Label supplied with bottle of concentrate
Graduated cylinders
Filtered Deionized Water
Anhydrous, Reagent Grade Isopropanol
Anhydrous, Reagent Grade Methanol
Anhydrous, Reagent Grade Acetone
Anhydrous, Reagent Grade Ethanol or Reagent Alcohol
Optional Supplies and/or Equipment
AC-038 – 5 L space-saver container

If bottles are reused be sure to clean the bottle thoroughly. Use 1/10 mixture of household bleach to disinfect the bottles.

WARNING! These reagents contain moderately toxic chemicals that require care in handling. Always use appropriate safety measures including gloves and eye protection when handling reagents. Refer to material safety data sheet (MSDS).

SS-141A Gram Reagent A, Safranin Concentrate prepared with Isopropanol and Acetone

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-041AA):

1. Obtain an empty 500 mL bottle.
2. Add 23 mL SS-141A to 382 mL Isopropanol and 95 mL Acetone.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-141A Gram Reagent A, Safranin Concentrate prepared with Isopropanol and Methanol

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-041A):

1. Obtain an empty 500 mL bottle.
2. Add 23 mL SS-141A to 287 mL Isopropanol and 190 mL Methanol.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-141AF Gram Reagent A, Fuchsin Concentrate prepared with Isopropanol and Acetone

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-041AAF):

1. Obtain an empty 500 mL bottle.
2. Add 23 mL SS-141AF to 382 mL Isopropanol and 95 mL Acetone.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-141AF Gram Reagent A, Fuchsin Concentrate prepared with Isopropanol and Methanol

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-041AAF):

1. Obtain an empty 500 mL bottle.
2. Add 23 mL SS-141AF to 287 mL Isopropanol and 190 mL Methanol.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-141B Gram Reagent B, Iodine Concentrate

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-041B):

1. Obtain an empty 500 mL bottle.
2. Add 50 mL SS-141B to 450 mL Deionized Water.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-141C Gram Reagent C, Crystal Violet Concentrate

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-041C):

1. Obtain an empty 500 mL bottle.
2. Add 15 mL SS-141C to 485 mL Deionized Water.
3. Mix well.
4. Label new bottle with supplied label.
5. Write LOT number and expiration date from concentrate bottle on new label.

SS-029C Nozzle Cleaning Solution Concentrate

Mixing Instructions (when mixed according to these instructions the product is equivalent to SS-029):

1. Add 250 mL Methanol or Ethanol to bottle of SS-029C.
2. Mix well.
3. Check box on label when mixed.

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