
NANODUCT[®] Patient Simulator AC-111



COPYRIGHT 2004 WESCOR INC

C 2004 Wescor, Inc. All rights reserved. Printed in the United States of America.

Wescor and Nanoduct are trademarks of Wescor, Inc.

Other trade names used in this manual may be trademarks of their respective owners, used here for information only.

U.S. Patent Number 6,198,953.

All information in this document is subject to change without prior notice.

section

1 INTRODUCTION

1.1	Customer Service
1.2	System Description

section

2 OPERATING THE NANODUCT PATIENT SIMULATOR

2.1	Setup
2.2	Operation
2.3	Cleanup
2.4	Battery Replacement

section

3 TROUBLESHOOTING

3.1	Troubleshooting	
-----	-----------------	--

APPENDIX A

Specifications	3	27
----------------	---	----

APPENDIX B

Supplies							
----------	--	--	--	--	--	--	--

i

1.1 Customer Service

Wescor is prepared to help you resolve any difficulty with the operation or performance of this device. If you are unable to solve a problem using the procedures described in this manual, please contact us.

Customers should contact Wescor by telephone, fax, or e-mail. Outside the U.S., many of our authorized dealers offer customer service and support.

TELEPHONE

435 752 6011

Extension 0 - Operator 171 - Orders

172 - Service

173 - Product Information and Pricing

TOLL FREE (US and Canada)

800 453 2725

Extension

- 0 Operator
- 171 Orders
- 172 Service
- 173 Product Information and Pricing

FAX 435 752 4127

E-MAIL service@wescor.com

WEBSITE:

www.wescor.com/biomedical

MAILING ADDRESS:

Wescor, Inc Biomedical Products Division 459 South Main Street Logan, Utah 84321-5294 USA

PLEASE NOTE:

Wescor's Authorized European Representative for matters relating to the Medical Device Directive is:



MT Promedt Consulting GmbH Altenhofstraße 80 D-66386 St. Ingbert, Germany Tel. +49 6894 581020 Fax +49 6894 581021 email info@mt-procons.com

1.2 System Description



Nanoduct Inducer/ Analyzer

Intended Use

The Nanoduct Patient Simulator is used as an accessory to the Nanoduct system by hospital or clinical personnel. The simulator allows measurements of prepared electrolyte solutions to be analyzed by the Nanoduct Neonatal Sweat Analysis System with the same degree of accuracy and in a similar manner to the total electrolyte concentration measurements made on the sweat from actual patients for the diagnosis of Cystic Fibrosis.



1.2 System Description



The Patient Simulator is a small DC-powered benchtop accessory with a receptacle where a standard Nanoduct Sensor can be located and locked into position. It is powered by a universal DC adapter, or can be powered by two standard 9 Volt batteries when Mains voltage is not available or convenient. An electronic circuit heats the sensor and holds it at approximately 31.5 degrees Celsius, simulating human body heat. When a liquid sample is introduced through the port on the front of the simulator, it enters the Nanoduct sensor through the apex of the bottom collector surface, just as a sweat sample from a patient would enter. The excess fluid from the sensor exit port runs into the waste collection pad and is easily cleaned using the provided tools.

The Patient Simulator can be used to measure unknown samples as part of a proficiency study or correlative experiments or to verify or validate the function and accuracy of the Nanoduct instrument. It is not intended to be a replacement for or alternative to the stable electronic calibration and control standards provided with the Nanoduct, which are less susceptible to degradation and human error, which can easily be verified to traceable standards, and which provide certainty about the function and accuracy of the Nanoduct electronic circuits. However, the simulator, used in conjunction with the calibration and control standards, allows independent confirmation of the Nanoduct's function, sample testing of the disposable Nanoduct sensors, and easy correlation with other analysis methods.



1.2 System Description



FRONT VIEW

1.2 System Description

INSTRUMENT TOP AND FRONT

Sensor Retaining Arm and Shroud

This is used to hold the Nanoduct sensor in place during simulation. The retaining arm is locked in place using the lever on the right side of the Simulator (see Right Side, below). The shroud helps minimize the effects of ambient air currents.

Sensor Receptacle

The raised round surface is designed to accept the Nanoduct sensor. The round hole at the bottom of the raised area is designed to accept waste collection pads to absorb excess fluid from the sensor output.

Ready Indicator

This LED indicator glows green when the Simulator has reached the acceptable temperature range (usually within 2 to 5 minutes of power-up when using the AC adapter. When using DC it depends of the state of the batteries).





TUBING FILL LENGTH

Power Indicator

This LED indicator glows yellow when power is applied or flashes yellow to indicate that batteries need to be replaced (battery mode only).

Restart Button

The instrument has an automatic shut-off at about 10 minutes to preserve battery life. Pressing the reset button renews this period and can be used to allow the instrument to reach optimum operating temperature inside the sensor cell.

"Tubing Fill Length" Scale

This graphically indicates the approximate length of tubing that should be filled with sample for optimum simulation results.



Tubing Input

The input nipple accepts the microbore tubing attached to the syringe for injecting sample solution into the Simulator. Note: Some Simulators have the rectangular shield shown at left.




Locking Lever







1.2 System Description



INSTRUMENT REAR

1.2 System Description



RIGHT SIDE Locking Lever

While pressing down on the top of the shroud with fingers on each side of the retainer button (avoid pressing on the retainer button), rotate the lever up and towards the back of the Simulator to lock the mounted sensor assembly into the Simulator.

INSTRUMENT REAR

ON/OFF Switch

- I = ON
- 0 = OFF

Universal DC Adapter Connector

Used to connect the Simulator to a universal DC power adapter.





Other Equipment and Supplies

Sensor (standard Nanoduct supply item) The base of the sensor is a shallow concavity leading at its center to an entry port, and from there to an internal fine channel passing by two analyzer micro-electrodes, forming a microconductivity cell.

Sensor Socket (part of the Nanoduct System)

The socket on the cable connects the electrode/sensor cable assembly to provide electrical current for sweat stimulation and to receive sensor signals to the analyzer module during the analysis phase.

Cleaning Adapter

The adapter is clamped into the Simulator in place of the sensor to allow the Simulator to be cleaned.

Note:

To preserve battery power when operating in battery mode, leave power off except while performing a test, (after the Average Conductivity is displayed and recorded). In addition, the Simulator shuts own automatically after approximately 10 minutes and must be turned on again to complete testing or retrieve a result.



1.2 System Description



Note:

The following are features of the Nanoduct Inducer Analyzer and not the Patient Simulator and are presented to allow effective use of the Simulator. Refer to the Nanoduct Neonatal Sweat Testing System User's Manual for more information.

DISPLAY (Nanoduct Inducer/Analyzer)

The LCD alphanumeric read-out displays all functions and results as they occur, including iontophoresis status and analysis results. The user responds to prompts on the display to activate various functions. The instrument emits a short beep once iontophoresis is complete.

The instrument displays concentration measurements in mmol/L (equivalent NaCl) during analysis. The concentration display is factory calibrated for a standard NaCl solution at 31.5 °C which has been established as the average temperature reached by the sensor microconductivity cell after 3 to 4 minutes in place on the patient's arm, calibrated to equivalent resistance.


LOW BATTERY Indicator (Nanoduct Inducer/Analyzer)

If the battery voltage drops below a preset level as the instrument is turned on, the Low Battery indicator appears on the display accompanied by a warning tone. When this indicator first sounds you can usually complete 1 to 2 tests before a dead battery is indicated. If the battery is too low to complete a test, the module will automatically shut off and iontophoresis cannot be started. You must replace the battery to continue. See Section 4.3 for information.

Note:

The inducer/analyzer is powered by one 9 Volt cell, which should provide a minimum of 150 tests with a standard alkaline battery or 300 tests with a lithium battery. SECTION 2 OPERATING THE PATIENT SIMULATOR

2.1 Setup




Inducer/Analyzer on a suitable lab bench, cart, or table. Place the units away from strong drafts in a room with normal ambient

If using the DC adapter, plug the adapter into the Simulator using the receptacle on the back of the unit. Plug the adapter into a

13

2.2 **Operation**



Snap a standard Nanoduct sensor into the sweat analysis connector of the Nanoduct with the exit hole on the same side as the

Attach a #22 blunt needle to the 1 mL syringe. Attach the provided length of microbore tubing to the blunt needle.

An ampule organizer (AC-071) is available for holding sweat control ampules. See Appendix B.

Open the sample to be measured and immerse the end of the microbore tubing in


Draw the sample fluid into the microbore tubing about 3.5 to 4.5 inches (about 90 to 115 mm). Sample volume should be about 35 microliters. The suggested tubing fill length is represented at the bottom of the front panel of the instrument to allow you to gauge the needed sample volume. If necessary, smaller volumes can be used.

Inspect the sample in the tubing for bubbles. Bubbles in the sample tubing will cause erratic readings. There must be at least 1 inch of bubble free sample near the free end of the tubing. If there is not a sufficient length of bubble-free sample, you must carefully repeat the process of drawing sample into the tubing until you achieve a sufficient length of bubble-free sample.


2.2 Operation



2.2 Operation



2.3 Cleanup



2.3 Cleanup



Attach the free end of the new tubing to the nipple of the cleaning adapter.

Place the cleaning adapter in the receptacle of the Simulator in the sensor position. Lock the retainer cup in place.

Push the deionized water through the system until the air in the tubing reaches the Simulator retaining cup. Work the water back and forth a few times with the plunger. Then push the plunger complete--ly down to expel all the remaining water with a burst of air.

Replace the 1 mL syringe with the 10 mL syringe and purge the Simulator with air. Fill the 10 mL syringe with air and push the air through the Simulator. Repeat at least 3 times. Remove the blunt needle from the 10 mL syringe after each purging in order to refill the 10 mL syringe with air.

Place a new waste collection pad in its holder below the sensor receptacle.

2.4 Battery Replacement

When used in battery mode, the Patient Simulator is powered by two 9 volt batteries. The batteries are mounted internally on a removable plate under the unit. To replace the batteries:



Completely loosen the two captive screws from the bottom plate.

Remove the bottom plate with the batteries attached.

B Remove both sp

Remove and discard (see following page) both spent batteries.



Insert 2 new 9 volt batteries (type NEDA/ANSI 1604A). Be sure the batteries are correctly aligned and seated in the connectors.

2.4 Battery Replacement



Battery Care and Disposal

WARNING!

Do not connect improperly, charge, or dispose in fire. Batteries may explode or leak. To prevent burn injury, do not carry batteries loose in your pocket.

Dispose of batteries in accordance with applicable laws and ordinances.

SECTION 3 TROUBLESHOOTING

SECTION 3 TROUBLESHOOTING

3.1 Troubleshooting

The Patient Simulator requires no routine maintenance except for cleaning as described in Section 2 and/or replacement of batteries when necessary. There are no user-serviceable parts.

SYMPTOMS

PROBLEM/SOLUTION

Temperature "ready"light does not come on.

AC adapter not properly plugged in. Dead batteries (without AC adapter). Possible calibration change. Contact your dealer or Wescor Service for calibration information.

Change in performance.

Check for stable ambient temperature: no cold drafts or excessive heat from other equipment.

Unexpected high reading when using control solutions.

Alternate conductivity path/damage to Simulator's sensor receptacle. Contact Wescor Service.

A P P E N D I X A SPECIFICATIONS

A P P E N D I X A SPECIFICATIONS

Power	7.5 to 14 Volts.
	Current less than .4 Amps during peak heating. Jack for external power supply: 12 V AC/DC adapter provided.
	Internal connections for two standard 9 Volt batteries (NEDA/ANSI 1604A). Auto shut off after 10 minutes. Can be extended by Restart push button.
Temperature	Holds liquid sample at 31.5 °C \pm 1 °C.
Keyboard	ON, OFF switch and RESTART button
Regulatory	Meets EN 61326 Standards for EMC compatibility. Meets IEC 60601 standards for safety. Manufactured under Wescor's quality system (ISO 13485). Bears the CE mark as a Class IIa medical device. Certificate issued by BSI 0086.
Size	Approximately 11.43 cm W x 10.92 cm D x 11.43 cm H. (4.5" x 4.3" x 4.5")
Weight	Approximately 0.7 kg (1.5 lbs.)

APPENDIX B SUPPLIES

A P P E N D I X B SUPPLIES

SUPPLIES

Nanoduct Sensor Cells	565-01
Ampule Organizer)71
DC Power Adapter AC-1	113
Cleaning Adapter	398-01
1 mL Tuberculin Syringe w/ #22 blunt needle SS-0)45
10 mL Tuberculin Syringe	175
Microbore Tubing lengths (14") SS-0)44
Cotton Swabs	290
Waste Collection PadsSS-2	220
Deionized Water SS-0)06
Hemostat Clamp	126
Plastic Tweezers	125