MYCOSCREEN <u>PLUS</u> is a methodiagnosis between commensal germs that could be incriminate infectious process (≥10 <sup>4</sup> CCU/ml) MYCOSCREEN <u>PLUS</u> must be a diagnosis based upon culture c	Non-gonoccocal urethritis +++- Urethro-prostatitis ++ Epididymitis ++ Urinary stones ++ Pyelonephritis - Pelvic inflammation - Post partum fever/abortion ++ Chorioamnionitis ++	dreary/ucurn and mycoprasma no encountered. U.u. and M.h. are sy pathogenic. Respiratory infections neonate as a result of contamin birth (3). In adults, the infections described in the table below (6):	capable of self-reproducing (0.1: species have been found to date class (derived from: <i>cutis mollis</i> ; s bacteria in their lack of a cell wall; to ß-lactams, as well as by the pr sterol obtained through their adhe Mycoplasmas are relatively fra acellular cultures in the presence at a constant temperature of 35 th Most human mycoplasmas are of that have been isolated from the	Mycoplasma hominis (M.h MYCOSCREEN <u>PLUS</u> must be us Agar plate for the enumerat identification of U.u. and M.h. Ear 2 - INTRODUCTION	1 - INTENDED USE MYCOSCREEN <u>PLUS</u> has been of differentiation of <i>Ureaplasma</i>	US-2007-09	MYCOSCR Detection and c of urogenital m 16 tests (Ret 32 tests (Ret
od that enables the differential germs (≤10 <sup>3</sup> CCU*/ml) and ed in an infectious or super t. *Colour Change Unit. ssociated with a conventional on A7 Agar plates, (modified	+ + + + + + + + + + + + + + + + + + + +	running are the most commonly exually transmitted and can be s or meningitis can occur in the ation from the genital tract at caused by U.u. and M.h. are realyticum M. hominis	5 to 0.25 $\mu$ m) (1). 15 humon all belonging to the mollicutae oft skin). They differ from other and hence a natural resistance resence of a membrane rich in ssion to eukaryotic cells. Since gile, they will only grow in of various growth factors and o 37°C. Commensal. Of the 9 species a urogenital tract, <i>Ureaplasma</i>	<ol> <li>in clinical specimens. sed in association with the A7 tion and the microscopical ch kit contains 16 or 32 tests.</li> </ol>	designed for the detection and a <i>urealyticum</i> (U.u.) and	(€	EEN <u>PLUS</u> Ifferentiation ycoplasmas f. 00024) f. 00025)
<ul> <li>6 - <u>SAMPLE COLLECTION</u></li> <li>6.1 <u>Sample collection</u></li> <li>6.1 <u>Sample collection</u></li> <li><u>Endocervical and vaginal sample collection</u></li> <li>Use only a Dacron® or rayon swab or a cytobrush to collect samples.</li> <li>Use only a bacron® or rayon swab or a cytobrush to collect samples.</li> <li>The cervix should be carefully cleaned with a swab, to remove secretions, before collecting the sample with a new swab. As</li> </ul>	<ul> <li>Reagents containing raw materials of animal origin must be handled with caution.</li> <li>Do not use reagents after the expiry date.</li> <li>Do not use reagents that have been damaged or that have been poorly conserved before use.</li> <li>A positive result with the MYCOSCREEN <u>PLUS</u> method indicates colonisation by urogenital mycoplasmas, but cannot alone be used to make a clinical diagnosis. This must be made by a doctor and is a function of the biological results and clinical signs.</li> </ul>	<ul> <li>PRECAUTIONS</li> <li>The reagents are intended solely for <i>in vitro</i> use and must be handled by authorized personnel.</li> <li>The patient samples and inoculated reagents are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.</li> </ul>	Reconstituted medium (UMMt+UMMIyo) -         Composition in g/l of distilled water:         Mycoplasma broth.       20         Foal serum.       200 ml         Yeast extract.       5,8         Cysteine.       0,3         Arginine.       3,6         Urea.       0,04         Phenol red.       0,04         Antibiotics.       10 ml         pH : 6.1 ± 0.1       0.1	S. M.h.: Mycoplama hominis growth activator (4.5ml)       1       1         UMMt medium - Composition in g/l of distilled water:       20         Mycoplasma broth	#00024       #00024       #00025         MYCOSCREEN <u>PLUS</u> tray: Divisible microplate       4       8         of 4 rows of 2 wells (U.u) / (M.h),       4       8         U.M.M.t: Vial of mycoplasma transport medium (2 ml)       16       32         U.M.M.lyo : Vial of lyophilised growth medium (+2ml UMMt)       16       32	4 - <u>REAGENTS</u> Description Quantity	<b>3 - PRINCIPLE</b> MYCOSCREEN <u>PLUS</u> is a 24-hour method based upon the ability of U.u. and M.h. to metabolise, respectively, urea and arginine. Mycoplasma growth in a liquid medium is visualised by the colour change of an indicator from yellow to fuchsia pink. This is due to the alkalinization of the medium as a result of the release of ammonia. MYCOSCREEN <u>PLUS</u> is carried out in a divisible microplate enabling the testing of 4 samples. Each microplate is composed of 4 rows of 2 wells: a <i>Ureaplasma urealyticum</i> (U.u) well containing lincomycin and a <i>Mycoplasma hominis</i> (M.h) well containing erythromycin.

Mycoplasmas adhere strongly to mucous cells, the mucous lining should be well scraped to obtain a rich specimen (2). Irethral sample collection

lean the meatus and swab or scrape the area to obtain cells.

ollect sperm or first micturition in a sterile tube or bottle. perm, Urine

### astric secretions

catheter and transferring to a sterile bottle. collect gastric secretions from the neonate by aspirating through a

## .2 Transport in UMMt medium

nomogenised liquid .iquid samples: Inoculate a vial of UMMt medium with 200µl of wab samples: Place the swab in a vial of UMMt medium.

Phe inoculated UMMt medium may be kept for 8 hours at room emperature (18-25°C) or 16 hours at 2-8°C.

# PREPARATION AND STORAGE OF REAGENTS

All the reagents are ready to use. The vials may be stored at 2-solution of the stored at 2- stored at the stored at 2- stored at 2

can be conserved in the original aluminium packet, sealed remetically, for up to 3 weeks at 2-8°C. vells are used, the remaining wells on the MYCOSCREEN <u>PLUS</u> tray When a single row of wells (U.u.) (M.h.) or two or three rows of

The Mh supplement is stable for 3 months after opening.

# - MATERIAL REQUIRED BUT NOT PROVIDED

- Sample collecting material (Swabs, cytobrushes, sterile
- containers for liquid samples)
- Pipettes and tips (100µl)
- Waste container for contaminated waste
- Paraffin oil
- A7 AGAR (#00079) or (#00090)

### - METHOD

### 1.1 Preparation of reagents

to be tested. Allow the reagents to reach room temperature Prepare the same number of rows of wells, as there are samples

he aid of the notches. If necessary, cut out one or several rows of wells (U.u)(M.h) with

Label each test.

§6.2). Seed the UMMt medium with the swab or 200 µl of liquid sample.

Transfer all of the UMMt medium to the UMMlyo vial. Do not lispose of the empty UMMt vial.
 Mix well

## 1.3 Inoculation of the tray wells

Distribute successively:

(U.u) well: (M.h) well: 100 µl of inoculated UMMlyo

50 µl of Mh supplement. 100 µl of inoculated UMMlyo

empty UMMt vial. Add two drops of paraffin oil to each of the wells. Remove 500  $\mu l$  of seeded UMMlyo medium and transfer it to the

ase of a positive screening/ detection result, an A7 Agar culture can e carried out (§ chapitre11). Store excess inoculated UMMIyo in its vial at 2-8°C, so that in the

Absence of mycoplasmas. <b>11 - QUALITY CONTROL</b> Quality control can be carried out from a lyophilised reference strain ( <i>Ureaplasma urealyticum</i> ATCC 27618). Prepare a preculture after inoculation of 100µl of the regenerated	Reincubate the tray and the vial. <b>10.2.2 Reading within 48 hours</b> <b>Yellow colour in the U.u and M.h wells and red colour in the</b> <b>UMMt vial</b> Probability that mycoplasmas at a rate of <10 <sup>2</sup> CCU/ml are present in the sample. <b>Yellow colour in the U.u and M.h wells and red colour in the</b> <b>UMMt vial</b>	(≥10 <sup>4</sup> CCU/ml). <u>Note</u> : Very strong concentrations of Uu or Mh can cause the medium in the two wells to change color. <b>Yellow colour in the U.u and M.h wells and red colour in the</b> <b>UMMt vial</b> Presence of mycoplasma at a rate of <10 <sup>4</sup> CCU/ml. In this case, seed an A7 Agar plate (§12) especially for samples other than the endocervical or urethral in origin (§13). <b>Yellow colour in the U.u and M.h wells and yellow colour in</b> <b>the UMMt vial</b> Absence of mycoplasmas, or presence of mycoplasmas at a rate of <10 <sup>5</sup> CCU/ml	An orangey colouring should be considered as a positive test. <b>10.2.1 Reading within 24 hours</b> <b>Red colour change in the U.u well:</b> Strong probability that <i>Ureaplasma urealyticum</i> at a supra- pathological rate is present in the sample ( $\geq$ 10 <sup>4</sup> CCU/ml). <b>Red colour change in the M.h well:</b> Strong probability that <i>Mycoplasma hominis</i> at a supra- pathological rate is present in the sample ( $\geq$ 10 <sup>4</sup> CCU/ml). <b>Red colour change in the U.u and M.h wells:</b> Strong probability that <i>Ureaplasma urealyticum</i> and <i>Mycoplasma</i> <i>hominis</i> at supra-pathological rates are present in the sample	<ul> <li>10.1 Validation Check that the 2 wells (U.u) (M.h) are limpid. A cloudy appearance indicates bacterial contamination. In this case, repeat the analysis.</li> <li>10.2 Reading and interpretation Urogenital Mycoplasma growth in the wells is indicated by the alkalinization of the medium, which changes colour from yellow to red.</li> <li>The medium remains yellow when no growth of urogenital Mycoplasma occurs.</li> </ul>	<ul> <li><u>Note</u>: Certain <i>M. hominis</i> strains exhibit a weak enzymatic activity. The addition of a drop of Mh Supplement allows a colour change to be visualized within 24 hours for strains whose concentration is ≥10<sup>4</sup>CCU/ml.</li> <li><b>9.4</b> Incubation Incubate tray wells and the UMMt vial at 37°C ± 1°C for 24 hours, then read the results. For enumeration, read the tray results within 24 hours. Prolong the incubation for up to 48 hours for the rates inferior to 10<sup>4</sup>CCU/ml or to enable strains with a weak enzymatic activity to become apparent.</li> <li><b>10 - READING AND INTERPRETATION</b></li> </ul>
<b>15 - PERFORMANCE</b> A study carried out on 334 clinical specimens (224 endocervical specimens, 25 urethral specimens, 75 sperm specimens, 5 urine specimens and 5 miscellaneous specimens), of which 99 specimens were positive, gave the following results compared with the	<ul> <li>change color.</li> <li>Some bacteria that are present in quantities of &gt;10<sup>6-7</sup> CFU/ml and contain urease may cause the two wells to change colour.</li> <li>A basic sample pH (pH &gt; 8) may lead the UMM medium to change colour.</li> <li>As for all germ search methods, the quality of the sample can influence the test result. A negative test does not therefore necessarily indicate the absence of infection.</li> </ul>	<ul> <li>13 - SPECIFIC CASES</li> <li>For sperm, urine, surgical salpingitis specimens, haemoculture, quality control of cell culture, etc, it is important to determine the concentration of mycoplasmas present. It is preferable that an A7 Agar plate is seeded with an inoculated UMMIyo medium, without waiting for the screening results.</li> <li>The concentration usually quoted for <i>U. urealyticum</i> in a first urine stream, sperm, or an endotracheal specimen is 10<sup>3</sup> CCU/ml (4).</li> <li>14 - LIMITATIONS OF THE PROCEDURE</li> <li>The presence of <i>Ureaplasma urealyticum</i> or <i>Mycoplasma hominis</i> at very high concentrations can cause the medium in the two wells to</li> </ul>	<ul> <li>12.2 Enumeration</li> <li>Enumeration is performed upon an upside down agar plate with the Enumeration is performed (x10 objective). The number of colonies per field of view is recorded:</li> <li>- Less than 1 coloniy / field: ≤10<sup>3</sup> CFU/ml</li> <li>- 1 to 5 colonies / field: approx. 10<sup>4</sup> CFU /ml</li> <li>- 5 to 10 colonies / field: approx. 10<sup>5</sup> CFU /ml</li> <li>- 10 to 20 colonies / field: approx. 10<sup>6</sup> CFU /ml</li> <li>- 20 colonies / field: approx. 10<sup>6</sup> CFU /ml</li> <li>CFU: Colony Forming Units</li> </ul>	<ul> <li>If the screening is positive, take the seeded UMMlyo medium conserved at 2-8°C and inoculate an A7 Agar plate with 100 µl of the homogenised suspension.</li> <li>Incubate at 37°C ± 1°C in an anaerobic environment for 48 hours.</li> <li>12.1 Morphological identification Ureaplasma urealyticum: Appearance of a brown-black precipitate (variable size, "sea urchin shaped"). The colonies are small.</li> <li>Mycoplasma hominis: "fried egg shaped" appearance. The colonies are larger than those of <i>U. urealyticum</i>.</li> </ul>	strain in a vial of reconstituted UMM medium (UMMt+UMMIyo). Incubate at $37^{\circ}C \pm 1^{\circ}C$ for 6 hours. Inoculate a new UMM broth reconstituted with 200µl of the homogenised preculture. Inoculate the two wells of the MYCOSCREEN <u>PLUS</u> tray and perform the test as indicated in these instructions (§9 et 10). <b>Expected results:</b> U.u M.h <b>A</b> Note: In the event of a lyophilised strain with low ureasic activity, incubation can be continued for up to 18-24 hours. The broth should be inoculated with 20 µl of preculture.

MYCOFAST *Evolution* 2 method, for rates of 10<sup>4</sup> CCU/ml. cellaneous specimens), of which 99 specimens on 334 clinical specimens (224 endocervical al specimens, 75 sperm specimens, 5 urine

otoction - 97 hours	II urosluticum	M hominic
Sensitivity	92.7%	100%
specificity	98.8%	99.3%
Pde	95.5%	91.6%
NPV	98.1%	100%
etection - 48 hours	U. urealyticum	M. hominis
Sensitivity	98.8%	95.7%
specificity	98.7%	98.2%
PA	96.6%	%00
UPV	99.5%	99.2%
Vith the MYCOSCREEN	<u>PLUS</u> tray 75,2%	of the U.urealy
traine and 18 30/ of the h	1 hominic strains w	ioro dotoctod off

2 tray. strains and 48,2% of the *M. nominis* strains were detected after 24 hours of incubation. For a detection threshold of 10<sup>4</sup>CCU/ml the results were similar to those obtained with the MYCOFAST *Evolution* ter 24 ticum

## Differentiation of U.u and M.h

MYCOFAST Evolution 2 method is 96.5%. between U. urealyticum and M. hominis, The global concordance for the differentiation within 24 hours compared to the

### 16 - WASTE ELIMINATION

and current regulations for this kind of product in the country of use Waste should be disposed of in accordance with the hygiene rules

### 17 - BIBLIOGRAPHY

Paris **1. BEBEAR C., DE BARBEYRAC B.** 1994. Les mycoplasmes, p. 1443-1463. Dans FRENEY J., RENAUD F. , HANSEN W. , BOLLET C. (éd.). "Manuel de bactériologie clinique", 2<sup>eme</sup> éd., **vol.3**, Elsevier,

génitaux. Feuil. Biol., 34: 21-24. BOUCAUD-MAITRE Y. et THOINET S. 1993. Analyse des prélèvements en bactériologie médicale - 2<sup>ème</sup> partie : prélèvements

molecular investigations of genital mycoplasmas from women and C., SEFFERT P. and POZZETTO B. 1995. Epidemiologic and 3. GRATTARD F., SOLEIHAC B., DE BARBEYRAC B., BEBEAR

4. PEREYRE S., BEBEAR C.M., BEBEAR C. 2001. Les

mycoplasmes en pathologie humaine. Revue Française Laboratoires. Supplément au N°**329**, 34-36. des

mycoplasmas) in primary cultures of clinical material. J. Clin. Microbiol. 3: 613-625. medium (A7) for identification of Ureaplasma urealyticum (human T 5. SHEPARD, M.C. and LUNCEFORD C.D.1976. Differential agar

Practices of Infectious Diseases, 4th ed., vol. 2, Churchill Mycoplasma) and *Mycoplasma hominis*, p. 1713-1718. Dans MANDELL G. L., BENNET J. E. and DOLIN R. (ed.). Principles and 6. TAYLOR-ROBINSON D. 1995. Ureaplasma urealyticum (T-strain Livingstone, New York.

Allée d'Athènes 83870 SIGNES FRANCE ELITech France Parc d'activités Place of manufacture :

