

FUNGICHROM® (#44325)

Tableau d'identification / Table of identification / Identifikationstabelle /
Tabella d'identificazione / Cuadro de identificación / Quadro de identificação

GENRE/GENUS	ESPÈCE, SPECIE, ART, ESPECIE, ESPÉCIE	GAL	PRO	ACT	ONPG	EPA	SGL	GLY	URE	POX	GAL-SAC	TRE	MAL	CEL	RAF	LAC
<i>Candida</i>	<i>albicans</i>	+	+	+	-	-	-	-	-	-	+	+	+	-	-	-
	<i>albicans</i> var. <i>stellatoidea</i>	+	-	+	-	-	-	-	-	-	+	V	+	-	-	-
	<i>lipolytica</i>	-	+	+	-	-	+	+	+	-	-	-	-	-	-	-
	<i>guilliermondii</i>	-	+	+ 80 %	-	-	-	V	-	-	+	+	+	+	+	-
	<i>lusitaniae</i>	-	+	-	-	-	-	-	-	-	+	+	+	+	-	-
	<i>parapsilosis</i>	-	+	-	-	-	-	+	-	-	(+)	(+)	(+)	-	-	-
	<i>inconspicua</i>	-	+	-	-	-	V	V	-	-	-	-	-	-	-	-
	<i>zeylanoides</i>	-	+	-	-	-	-	-	-	-	-	V	-	-	-	-
	<i>famata</i>	-	V	-	-	-	-	-	V	-	+	+	+	+	+	-
	<i>kefyr</i>	-	-	+	+	-	-	-	-	-	+	-	-	-	+	+
	<i>krusei</i>	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
	<i>glabrata</i>	-	-	-	-	+ 97 %	V	V	-	-	-	+	-	-	-	-
	<i>rugosa</i>	-	-	-	-	-	+	-	-	-	+	-	-	-	-	-
	<i>tropicalis</i>	+ 8 %	-	V	-	-	V	+	-	-	+	+	+	-	-	-
<i>Saccharomyces</i>	<i>cerevisiae</i>	-	-	V	-	-	-	-	-	-	+	V	+	-	+	-
<i>Rhodotorula</i>	<i>rubra</i>	-	+	-	-	-	-	-	+	-	+	V	V	V	+	-
<i>Cryptococcus</i>	<i>neoformans</i>	-	-	-	-	-	-	-	+	+	+	V	+	V	V	-
	<i>uniguttulatus</i>	-	+	-	-	-	-	-	+	-	+	+	+	-	-	-
	<i>albidus</i>	-	+	-	-	-	-	-	+	-	+	+	+	+	V	-
	<i>laurentii</i>	-	-	-	-	-	+	+	+	-	+	+	+	+	+	+
<i>Trichosporon</i>	<i>cutaneum</i>	V	-	-	+ 60 %	-	V	V	+ 70 %	-	+	+	+	+	V	V
<i>Geotrichum</i>	<i>candidum</i>	V	-	+	-	V	V	V	-	-	+ 80 %	-	-	-	-	-
	<i>capitatum</i>	-	+	+	-	V	+	+	V	-	-	-	-	-	-	-
<i>Kloeckera</i>	<i>apis / apiculata</i>	-	-	+	-	-	-	-	-	-	-	-	-	+	-	-

+ : Caractère positif / Positive characteristic / Positiver Charakter / Carattere positivo / Carácter positivo

(+) : Caractère généralement positif / In general, positive characteristic / Im allgemeinen positiver Charakter / In generale, carattere positivo / En general, carácter positivo / Em geral, carácter positivo

- : Caractère négatif / Negative characteristic / Negativer Charakter / Carattere negativo / Carácter negativo

■ : Caractère majeur / predominant characteristic / Hauptcharakter / carattere principale / carácter principal / carácter essencial

FUNGICHROM® (#44325)

Tableau d'identification par code / Table of identification by code / Identifikationstabelle durch Code /
Tabella d'identificazione per codice / Cuadro de identificación por código / Quadro de identificação por código

CODE, CODICE, CÓDIGO	ESPÈCE, SPECIE, ART, ESPECIE, ESPÉCIE	CODE, CODICE, CÓDIGO	ESPÈCE, SPECIE, ART, ESPECIE, ESPÉCIE	CODE, CODICE, CÓDIGO	ESPÈCE, SPECIE, ART, ESPECIE, ESPÉCIE
00 00 20	<i>Candida glabrata</i>	20 00 00*	<i>Candida inconspicua</i>	40 00 01	<i>Kloeckera apis/apiculata</i>
00 00 52	<i>Saccharomyces cerevisiae</i>	20 00 00*	<i>Candida zeylanoides</i>	40 00 10	<i>Geotrichum candidum</i>
00 00 72	<i>Saccharomyces cerevisiae</i>	20 00 01	<i>Candida inconspicua</i>	40 00 52	<i>Saccharomyces cerevisiae</i>
00 01 51	<i>Trichosporon cutaneum</i>	20 00 02	<i>Candida guilliermondii</i>	40 00 72	<i>Saccharomyces cerevisiae</i>
00 01 73*	<i>Trichosporon cutaneum</i>	20 00 41	<i>Candida lusitaniae</i>	40 20 70	<i>Candida tropicalis</i>
00 01 73*	<i>Candida famata</i>	20 00 61	<i>Candida lusitaniae</i>	40 30 70	<i>Candida tropicalis</i>
00 03 13	<i>Cryptococcus neoformans</i>	20 00 70	<i>Candida parapsilosis</i>	41 00 10	<i>Geotrichum candidum</i>
00 03 43	<i>Cryptococcus neoformans</i>	20 00 71	<i>Candida lusitaniae</i>	41 00 16	<i>Candida kefir</i>
00 03 52	<i>Cryptococcus neoformans</i>	20 00 73*	<i>Candida guilliermondii</i>	41 10 16	<i>Candida kefir</i>
00 03 62	<i>Cryptococcus neoformans</i>	20 00 73*	<i>Candida famata</i>	41 10 56	<i>Candida kefir</i>
00 03 73	<i>Cryptococcus neoformans</i>	20 01 52	<i>Rhodotorula rubra</i>	41 30 10	<i>Geotrichum candidum</i>
00 10 00*	<i>Candida krusei</i>	20 01 70	<i>Cryptococcus uniguttulatus</i>	42 00 00	<i>Geotrichum candidum</i>
00 10 00*	<i>Candida rugosa</i>	20 01 72	<i>Rhodotorula rubra</i>	42 00 10	<i>Geotrichum candidum</i>
00 10 10*	<i>Candida krusei</i>	20 01 73*	<i>Candida famata</i>	42 30 00	<i>Geotrichum candidum</i>
00 10 10*	<i>Candida rugosa</i>	20 01 73*	<i>Rhodotorula rubra</i>	42 30 10	<i>Geotrichum candidum</i>
00 10 12	<i>Candida krusei</i>	20 01 73*	<i>Cryptococcus albidus</i>		
00 11 55	<i>Trichosporon cutaneum</i>			50 20 50	<i>Candida stellatoidea</i>
00 20 70	<i>Candida tropicalis</i>	20 10 00	<i>Candida inconspicua</i>	52 00 10	<i>Geotrichum candidum</i>
00 30 20	<i>Candida glabrata</i>	20 20 00	<i>Candida parapsilosis</i>	52 30 10	<i>Geotrichum candidum</i>
00 30 70	<i>Candida tropicalis</i>	20 20 12	<i>Candida guilliermondii</i>		
00 31 65	<i>Trichosporon cutaneum</i>	20 20 20	<i>Candida parapsilosis</i>	60 00 73	<i>Candida guilliermondii</i>
00 31 75	<i>Trichosporon cutaneum</i>	20 20 32	<i>Candida guilliermondii</i>	60 20 73*	<i>Candida guilliermondii</i>
00 31 77	<i>Cryptococcus laurentii</i>	20 20 33	<i>Candida guilliermondii</i>	60 20 73*	<i>Candida famata</i>
01 00 73	<i>Trichosporon cutaneum</i>	20 20 40	<i>Candida parapsilosis</i>	60 30 73	<i>Candida guilliermondii</i>
01 01 73	<i>Trichosporon cutaneum</i>	20 20 50	<i>Candida parapsilosis</i>	60 31 00*	<i>Geotrichum capitatum</i>
01 31 71	<i>Trichosporon cutaneum</i>	20 20 60	<i>Candida parapsilosis</i>	60 31 00*	<i>Candida lipolytica</i>
01 31 77	<i>Trichosporon cutaneum</i>	20 20 70	<i>Candida parapsilosis</i>	62 30 00	<i>Geotrichum capitatum</i>
01 33 77	<i>Trichosporon cutaneum</i>	20 20 73	<i>Candida guilliermondii</i>		
01 60 77	<i>Trichosporon cutaneum</i>	20 30 00	<i>Candida inconspicua</i>	70 00 50	<i>Candida albicans</i>
02 00 20	<i>Candida glabrata</i>	20 31 12	<i>Rhodotorula rubra</i>	70 00 60	<i>Candida albicans</i>
02 00 72	<i>Candida glabrata</i>	20 31 52	<i>Rhodotorula rubra</i>	70 00 70	<i>Candida albicans</i>
02 10 20	<i>Candida glabrata</i>	20 31 75	<i>Trichosporon cutaneum</i>	70 20 70	<i>Candida albicans</i>
02 30 20	<i>Candida glabrata</i>				
10 20 70	<i>Candida tropicalis</i>				
10 20 72	<i>Candida tropicalis</i>				
11 01 00	<i>Trichosporon cutaneum</i>				
11 01 71	<i>Trichosporon cutaneum</i>				
11 01 73	<i>Trichosporon cutaneum</i>				
11 01 77	<i>Trichosporon cutaneum</i>				
11 30 71	<i>Trichosporon cutaneum</i>				
11 30 72	<i>Trichosporon cutaneum</i>				
11 31 31	<i>Trichosporon cutaneum</i>				
11 31 45	<i>Trichosporon cutaneum</i>				
11 31 55	<i>Trichosporon cutaneum</i>				
11 31 71	<i>Trichosporon cutaneum</i>				