

ACID FAST REAGENT CKY	SÄUREFEST REAGENZ CKY
Carbol Fuchsin (Kinyoun) SS-61CKY	Karbofuchsin (Kinyoun) SS-061CKY
For in vitro diagnostic use with the Aerospray® Slide Stainer for acid fast bacteria staining.	In Vitro Diagnostikum für denAerospray® Färbeautomat zum Färben säurefester Bakterien.
WARNING: TOXIC, SUSPECTED CARCINOGEN. May be harmful if swallowed. Use only with adequate ventilation. Wear chemical protection for eyes and skin. Store at room temperature away from direct sunlight. Do not puncture package-contents will stain.	WARNUNG: GIFTIG, STEHT IM VERDACHT KREBS ZU ERREGEN. Verschlucken kann gesundheitsschädlich sein. Nur bei ausreichender Lüftung verwenden. Chemische Schutzkleidung für Augen und Haut tragen. Bei Zimmertemperatur ohne direkte Sonneneinstrahlung lagern. Packung nicht beschädigen – Inhalt färbt.
CONTENTS (500 mL): deionized water, ethanol, phenol, and basic fuchsin.	INHALTSSTOFFE (500 ML): entionisiertes Wasser, Ethanol, Phenol und basisches Fuchsin.
LOT EXP	CHARGE VERFALLDATUM