

ACID FAST REAGENT BP	SÄUREFEST REAGENZ BP
Potassium Permanganate SS-061BP	Kaliumpermanganat SS-061BP
For in vitro diagnostic use with the Aerospray® Slide Stainer for fluorescent staining of acid fast bacteria.	In Vitro Diagnostikum für den Aerospray® Färbeautomat zum Fluoreszenzfärben säurefester Bakterien.
WARNING: TOXIC, IRRITANT, OXIDANT, DO NOT FREEZE. Store at room temperature away from direct sunlight. Do not puncture package-contents will stain.	WARNUNG: GIFTIG, REIZEND, NICHT EINFRIEREN. Bei Zimmertemperatur ohne direkte Sonneneinstrahlung lagern. Packung nicht beschädigen – Inhalt färbt.
CONTENTS (500 mL): deionized water and potassium permanganate.	INHALTSSTOFFE (500 ML): entionisiertes Wasser und Kaliumpermanganat.
LOT EXP	CHARGE VERFALLDATUM