

HEMATOLOGY REAGENT A	HÄMATOLOGIE REAGENZ A
Basofix™ Rinse SS-149A	Basofix™ Spülung SS-149A
For in vitro diagnostic use with the Aerospray® Slide Stainer for hematology staining.	In Vitro Diagnostikum für den Aerospray® Färbeautomat zum Färben hämatologischen Materials.
WARNING: TOXIC, DO NOT FREEZE Vapor harmful; avoid prolonged or repeated breathing of vapor. Harmful if swallowed. Causes eye irritation. Store at room temperature away from direct sunlight. Do not puncture package-contents will stain.	WARNUNG: GIFTIG, NICHT EINFRIEREN. Dampf gesundheitsschädlich; längeres oder wiederholtes Einatmen vermeiden. Verschlucken ist gesundheitsschädlich. Reizt die Augen. Bei Zimmertemperatur ohne direkte Sonneneinstrahlung lagern. Packung nicht beschädigen – Inhalt färbt.
CONTENTS (if properly diluted): deionized water, SD3A ethanol, methanol, phosphate buffer, eosin dye, and preservative.	INHALTSSTOFFE (bei korrekter Verdünnung): entionisiertes Wasser, SD3A Ethanol, Methanol, Phosphatpuffer, Eosinfarbstoff und Konservierungsmittel.
LOT EXP	CHARGE VERFALLDATUM
PREPARED BY DATE	AUFBEREITET VON DATUM