


Para la calibración del analizador[™]
de conductividad del sudor Sweat-Chek[™]
Contiene: 60 ampollas conteniendo 0,75 mL
por ampolla de solución ClNa. Molaridad 90
mmol/L



For calibrating the Sweat-Chek[™]
Sweat Conductivity Analyzer
Contents: 60 ampules containing 0.75 mL
per ampule of NaCl solution. Molarity 90 mmol/L.
For in vitro diagnostic use.
Manufactured for Wescor, Inc.

WESCOR[®]



For calibrating the Sweat-Chek[™] Sweat Conductivity Analyzer

Contents: 60 ampules containing 0.75 mL per ampule of NaCl solution,
Molarity 90 mmol/L

For in vitro diagnostic use.
Manufactured for WESCOR, INC
Order Number SS-140

©1996 WESCOR, INC 450 South Main Street, Logan, Utah 84321 Made in the USA

Para la calibración del analizador de conductividad
del sudor Sweat-Chek[™]

Contiene: 60 ampollas conteniendo 0,75 mL por
ampolla de solución ClNa. Molaridad 90 mmol/L

Para uso en diagnóstico in vitro
Fabricado por WESCOR, INC.

Número de pedido SS-140



For calibrating the Sweat-Chek[™]

Sweat

Conductivity Analyzer

Contents: 60 ampules containing 0.75 mL
per ampule of NaCl solution. Molarity 90 mmol/L.
For in vitro diagnostic use.
Manufactured for Wescor, Inc.

Reorder Number SS-140

Exp. Date

Lot Number

90 mmol/L

Reorder Number SS-140

Exp. Date

Lot Number

© 1996 WESCOR, INC
450 South Main Street, Logan, Utah 84321

Made in the USA

L79A75

Nº de pedido SS-140
Fecha de caducidad
Nº de lote