

GRAM REAGENT C	GRAM REAGENZ C
Crystal Violet SS-041C	Kristallviolett SS-041C
For in vitro diagnostic use with the Aerospray® Slide Stainer for Gram staining.	In Vitro Diagnostikum für den Aerospray® Färbeautomat zum Gram- Färben.
WARNING: TOXIC, IRRITANT, POSSIBLE MUTAGEN, SUSPECTED CARCINOGEN, DO NOT FREEZE Store at room temperature away from direct sunlight. Do not puncture package – contents will stain.	WARNUNG: GIFTIG, REIZSTOFF, MÖGLICHES MUTAGEN, STEHT IM VERDACHT KREBS ZU ERREGEN, NICHT EINFRIEREN Bei Zimmertemperatur ohne direkte Sonneneinstrahlung lagern. Packung nicht beschädigen – Inhalt färbt.
CONTENTS (500 mL): deionized water, crystal violet, preservative, and surfactant.	INHALT (500 ML): entionisiertes Wasser, Kristallviolett, Konservierungsmittel und Tensid.