

For calibrating the Sweat-Chek™ Conductivity Analyzer	Zur Kalibrierung des Sweat-Chek™ Leitfähigkeitsanalysator
Contents: 60 ampules containing 0.75 mL per ampule of NaCl solution, Molarity 90 mmol/L.	Inhalt: 60 Ampullen mit je 0,75 ml Natriumchloridlösung, Molarität 90 mmol/L.
For in vitro diagnostic use.	In Vitro Diagnostikum.
Manufactured for WESCOR, INC	Hergestellt für WESCOR, INC