



1000 mmol/kg OPTI-MOLE™ **WESCOR**
 For calibration of Vapor Pressure Osmometers or Freezing Point Osmometers.

OSMOLALITY STANDARD
 Contents: 60 ampules containing 0.4 mL sodium chloride solution.
 For in vitro diagnostic use.
 Manufactured for Wescor, Inc.

PRINT-0104-01A

	For Calibration of vapor pressure osmometers or freezing point osmometers.	Zur Kalibrierung von Dampfdruckosmometern und Gefrierpunktosmometern.
	Osmolality standard	Osmolalitätsstandard
	Contents: 60 ampules containing 0.4 mL sodium chloride solution.	Inhalt: 60 Ampullen mit 0,4 ml Natriumchloridlösung.
	For in vitro diagnostic use.	In Vitro Diagnostikum.
	Manufactured for Wescor, Inc.	Hergestellt für Wescor, Inc.